

## Silver nanoparticle SERS substrate Instructions for Use

Nikalyte silver nanoparticle SERS substrates are supplied in packs of 5 paper substrates packed in a glass vial purged with Argon. The substrate can either be used [as supplied](#) (see instructions for use of unmounted SERS substrate on p2) or [mounted on the glass slides](#), which are also provided in the pack (see instructions for mounting and using the SERS substrates on glass slides on p3)





Figure 1 Nikalyte Silver SERS substrates kit

<b>Specification</b>	
<b>SERS active material</b>	Silver nanoparticles
<b>Sensitivity</b>	ppm to ppb
<b>Laser Wavelength</b>	532nm, 785nm
<b>Dimensions</b>	75mm x 25mm
<b>Active area</b>	8mm x 8mm
<b>Analyte application area</b>	4mm x 4mm
<b>Lifetime (unopened)</b>	3 months
<b>Lifetime (open pack)</b>	1 week

## Procedure for as supplied SERS substrates

- Remove paper substrate from the vial without touching the active area to avoid contamination.
- Do not wash the substrate or expose to ultrasound.
- Place paper substrate onto a clean absorbent tissue, such as filter paper or low lint paper towel, to absorb any excess analyte.
- Slowly and carefully inject analyte solution directly onto the uncoated area of the paper.
- Bring pipette into contact with the bare substrate close to the coated area.
- Apply sufficient analyte to just wet the substrate active area but avoid having non-absorbed liquid on the substrate surface as excess analyte liquid can reduce the Raman signal. 20 to 30  $\mu\text{l}$  of analyte is usually sufficient.
- Avoid dropping or splashing the analyte from above.
- Do not dip the substrate into the analyte.
- Measure the Raman signal immediately.
- Recommended laser power density  $15\text{W}/\text{cm}^2$
- Maximum laser power density  $20\text{W}/\text{cm}^2$



## Procedure for glass mounted SERS substrates

- Peel off the protective layer from the carbon sticky pad on the glass slide.
- Remove paper substrate from the vial with a clean pair of tweezers without touching the active area to avoid contamination.
- Do not wash the substrate or expose to ultrasound.
- Place the paper substrate onto the carbon pad and press gently on the corners to ensure the pad is stuck securely.
- Slowly and carefully inject analyte solution directly on to the uncoated area of the paper.
- Bring pipette into contact with the bare substrate close to the coated area.
- Apply sufficient analyte to just wet the substrate active area but avoid having non-absorbed liquid on the substrate surface as excess analyte liquid can reduce the Raman signal. 20 to 25  $\mu\text{l}$  of analyte is normally sufficient.
- Avoid dropping or splashing the analyte from above.
- Do not dip the substrate into the analyte.
- Measure the Raman signal immediately.
- Recommended laser power density  $15\text{W}/\text{cm}^2$
- Maximum laser power density  $20\text{W}/\text{cm}^2$

